

Anti-NRG4 antibody (11-60) (STJ192544)

STJ192544

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Pro-Neuregulin-4-Membrane-Bound Isoform (11-60) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	145957
Gene Symbol	NRG4
Uniprot ID	NRG4_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 11-60
Immunogen Region	11-60
Specificity	NRG4 polyclonal antibody (Pro-Neuregulin-4-Membrane-Bound Isoform) binds to endogenous Pro-Neuregulin-4-Membrane-Bound Isoform at the amino acid region 11-60.
Immunogen Sequence	



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4°C over night

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081