

Anti-PZP antibody (1301-1350) (STJ192551)

STJ192551

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Pregnancy Zone Protein (1301-1350) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	PZP
Uniprot ID	PZP_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 1301-1350
Immunogen Region	1301-1350
Specificity	PZP polyclonal antibody (Pregnancy Zone Protein) binds to endogenous Pregnancy Zone Protein at the amino acid region 1301-1350.
Immunogen Sequence	