

Anti-TGFBRAP1 antibody (671-720) (STJ192560)

STJ192560

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Transforming Growth Factor-Beta Receptor-Associated Protein 1 (671-720) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9392
Gene Symbol	TGFBRAP1
Uniprot ID	TGFA1_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 671-720
Immunogen Region	671-720
Specificity	TGFBRAP1 polyclonal antibody (Transforming Growth Factor-Beta Receptor-Associated Protein 1) binds to endogenous Transforming Growth Factor-Beta Receptor-Associated Protein 1 at the amino acid region 671-720.
Immunogen Sequence	