

## Anti-SPG7 antibody (71-120) (STJ192564)

STJ192564

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Paraplegin (71-120) is suitable for use in Western Blot and ELISA research applications.
<b>Applications</b>	WB, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Formulation</b>	PBS, 50% Glycerol and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	6687
<b>Gene Symbol</b>	SPG7
<b>Uniprot ID</b>	SPG7_HUMAN
<b>Immunogen</b>	Synthesized peptide derived from human protein at aa range 71-120
<b>Immunogen Region</b>	71-120
<b>Specificity</b>	SPG7 polyclonal antibody (Paraplegin) binds to endogenous Paraplegin at the amino acid region 71-120.
<b>Immunogen Sequence</b>	