

Anti-LDB3 antibody (41-90) (STJ192577)

STJ192577

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Lim Domain-Binding Protein 3 (41-90) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11155
Gene Symbol	LDB3
Uniprot ID	LDB3_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 41-90
Immunogen Region	41-90
Specificity	LDB3 polyclonal antibody (Lim Domain-Binding Protein 3) binds to endogenous Lim Domain-Binding Protein 3 at the amino acid region 41-90.
Immunogen Sequence	