

Anti-GIGYF2 antibody (1071-1120) (STJ192597)

STJ192597

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Grb10-Interacting Gyf Protein 2 (1071-1120) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26058
Gene Symbol	GIGYF2
Uniprot ID	GGYF2_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 1071-1120
Immunogen Region	1071-1120
Specificity	GIGYF2 polyclonal antibody (Grb10-Interacting Gyf Protein 2) binds to endogenous Grb10-Interacting Gyf Protein 2 at the amino acid region 1071-1120.
Immunogen Sequence	