

Anti-SACS antibody (4291-4340) (STJ192601)

STJ192601

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Sacsin (4291-4340) is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-300
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26278
Gene Symbol	SACS
Uniprot ID	SACS_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 4291-4340
Immunogen Region	4291-4340
Specificity	SACS polyclonal antibody (Sacsin) binds to endogenous Sacsin at the amino acid region 4291-4340.
Immunogen Sequence	