

Anti-ARAP1 antibody (STJ192917)

STJ192917

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Arf-Gap With Rho-Gap Domain-Ank Repeat And Ph Domain-Containing Protein 1 is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116985
Gene Symbol	ARAP1
Uniprot ID	ARAP1_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Immunogen Region	
Specificity	ARAP1 polyclonal antibody (Arf-Gap With Rho-Gap Domain-Ank Repeat And Ph Domain-Containing Protein 1) binds to endogenous Arf-Gap With Rho-Gap Domain-Ank Repeat And Ph Domain-Containing Protein 1.
Immunogen Sequence	