

Anti-ANKS1A antibody (STJ192929)

STJ192929

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ankyrin Repeat And Sam Domain-Containing Protein 1a is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23294
Gene Symbol	ANKS1A
Uniprot ID	ANS1A_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Immunogen Region	
Specificity	ANKS1A polyclonal antibody (Ankyrin Repeat And Sam Domain-Containing Protein 1a) binds to endogenous Ankyrin Repeat And Sam Domain-Containing Protein 1a.
Immunogen Sequence	