

Anti-DYNC1LI1 antibody (STJ192934)

STJ192934

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cytoplasmic Dynein 1 Light Intermediate Chain 1 is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	51143
Gene Symbol	DYNC1LI1
Uniprot ID	DC1L1_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Immunogen Region	
Specificity	DYNC1LI1 polyclonal antibody (Cytoplasmic Dynein 1 Light Intermediate Chain 1) binds to endogenous Cytoplasmic Dynein 1 Light Intermediate Chain 1.
Immunogen Sequence	