

Anti-PHLPP1 antibody (STJ193083)

STJ193083

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ph Domain Leucine-Rich Repeat-Containing Protein Phosphatase 1 is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23239
Gene Symbol	PHLPP1
Uniprot ID	PHLPP1_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Immunogen Region	
Specificity	PHLPP1 polyclonal antibody (Ph Domain Leucine-Rich Repeat-Containing Protein Phosphatase 1) binds to endogenous Ph Domain Leucine-Rich Repeat-Containing Protein Phosphatase 1.
Immunogen Sequence	