

Anti-DIAPH2 antibody (STJ193155)

STJ193155

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein Diaphanous Homolog 2 is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-P, IF-P, WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1730
Gene Symbol	DIAPH2
Uniprot ID	DIAP2_HUMAN
Immunogen	Synthesized peptide derived from human DIAP2
Immunogen Region	
Specificity	DIAPH2 polyclonal antibody (Protein Diaphanous Homolog 2) binds to endogenous Protein Diaphanous Homolog 2.
Immunogen Sequence	