

Anti-PRCP antibody (STJ193169)

STJ193169

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Lysosomal Pro-X Carboxypeptidase is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-P, IF-P, WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5547
Gene Symbol	PRCP
Uniprot ID	PCP_HUMAN
Immunogen	Synthesized peptide derived from human PCP
Immunogen Region	
Specificity	PRCP polyclonal antibody (Lysosomal Pro-X Carboxypeptidase) binds to endogenous Lysosomal Pro-X Carboxypeptidase.
Immunogen Sequence	