

Anti-XRCC5 antibody (STJ193184)

STJ193184

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-X-Ray Repair Cross-Complementing Protein 5 is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-P, IF-P, WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7520
Gene Symbol	XRCC5
Uniprot ID	XRCC5_HUMAN
Immunogen	Synthesized peptide derived from human Ku-80
Immunogen Region	
Specificity	XRCC5 polyclonal antibody (X-Ray Repair Cross-Complementing Protein 5) binds to endogenous X-Ray Repair Cross-Complementing Protein 5.
Immunogen Sequence	