

## Anti-FLT3 antibody (STJ193204)

STJ193204

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Receptor-Type Tyrosine-Protein Kinase Flt3 is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications. |
| <b>Applications</b>      | IHC-P, IF-P, WB  |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human, Rat, Mouse  |

### PRODUCT PROPERTIES

|                       |  |
|-----------------------|--|
| <b>Clonality</b>      | Polyclonal   |
| <b>Clone ID</b>       |  |
| <b>Concentration</b>  | 1 mg/mL  |
| <b>Conjugation</b>    | Unconjugated   |
| <b>Purification</b>   | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.          |
| <b>Dilution Range</b> | IHC-P 1:50-200<br>WB 1:500-2000  |
| <b>Formulation</b>    | PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.  |
| <b>Isotype</b>        | IgG  |
| <b>Storage</b>        | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |
| <b>Instruction</b>    |  |

### TARGET INFORMATION

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | <a href="#">2322</a>  |
| <b>Gene Symbol</b>        | <a href="#">FLT3</a>  |
| <b>Uniprot ID</b>         | <a href="#">FLT3_HUMAN</a>  |
| <b>Immunogen</b>          | Synthesized peptide derived from human Flt3   |
| <b>Immunogen Region</b>   |   |
| <b>Specificity</b>        | FLT3 polyclonal antibody (Receptor-Type Tyrosine-Protein Kinase Flt3) binds to endogenous Receptor-Type Tyrosine-Protein Kinase Flt3. |
| <b>Immunogen Sequence</b> |   |