

Anti-Phospho-NTRK2-Y817 antibody (STJ193227)

STJ193227

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Bdnf/Nt-3 Growth Factors Receptor-Y817 is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-P, IF-P, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4915
Gene Symbol	NTRK2
Uniprot ID	NTRK2_HUMAN
Immunogen	Synthesized peptide derived from human Trk B (Phospho-Y817) Polyclonal Antibody
Immunogen Region	
Specificity	Phospho-NTRK2-Y817 polyclonal antibody (Bdnf/Nt-3 Growth Factors Receptor) binds to endogenous Bdnf/Nt-3 Growth Factors Receptor only when phosphorylated at Y817.
Immunogen Sequence	