

Anti-SARS-nCoV-2 antibody [S03FC] (STJ193234)

STJ193234

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-SARS-nCov-2 Spike is suitable for use in ELISA, ECL and Immunogold research applications.
Applications	ELISA, ECL, Immunogold
Host/Source	Human
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S03FC
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	ELISA 1:5000-1:100000 WB 1:1000-1:10000
Formulation	PBS
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	The recombinant SARS-CoV-2 (2019-nCoV) Spike Protein S
Immunogen Region	
Specificity	This antibody detects SARS-CoV-2 (2019-nCoV) Spike Protein S1 domain
Immunogen Sequence	