

Anti-COLGALT1 antibody (STJ193464)

STJ193464

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Procollagen Galactosyltransferase 1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79709
Gene Symbol	COLGALT1
Uniprot ID	GT251_HUMAN
Immunogen	Synthesized peptide derived from human COLGALT1
Immunogen Region	
Specificity	COLGALT1 polyclonal antibody (Procollagen Galactosyltransferase 1) binds to endogenous Procollagen Galactosyltransferase 1.
Immunogen Sequence	



Western blot analysis of lysates from PC-12 cells, primary antibody was diluted at 1:1000, 4°C over night