

Anti-UAP1L1 antibody (STJ193968)

STJ193968

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Udp-N-Acetylhexosamine Pyrophosphorylase-Like Protein 1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	91373
Gene Symbol	UAP1L1
Uniprot ID	UAP1L_HUMAN
Immunogen	Synthesized peptide derived from human UAP1L
Immunogen Region	
Specificity	UAP1L1 polyclonal antibody (Udp-N-Acetylhexosamine Pyrophosphorylase-Like Protein 1) binds to endogenous Udp-N-Acetylhexosamine Pyrophosphorylase-Like Protein 1.
Immunogen Sequence	



Western blot analysis of lysates from AD293 cells. primary antibody was diluted at 1:1000, 4°C over night

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081