

Anti-TTC29 antibody (STJ194557)

STJ194557

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tetratricopeptide Repeat Protein 29 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	83894
Gene Symbol	TTC29
Uniprot ID	TTC29_HUMAN
Immunogen	Synthesized peptide derived from human TTC29
Immunogen Region	
Specificity	TTC29 polyclonal antibody (Tetratricopeptide Repeat Protein 29) binds to endogenous Tetratricopeptide Repeat Protein 29.
Immunogen Sequence	



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4°C over night