

Anti-GALNTL5 antibody (STJ194558)

STJ194558

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Inactive Polypeptide N-Acetylgalactosaminyltransferase-Like Protein 5 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	168391
Gene Symbol	GALNTL5
Uniprot ID	GLTL5_HUMAN
Immunogen	Synthesized peptide derived from human GLTL5
Immunogen Region	
Specificity	GALNTL5 polyclonal antibody (Inactive Polypeptide N-Acetylgalactosaminyltransferase-Like Protein 5) binds to endogenous Inactive Polypeptide N-Acetylgalactosaminyltransferase-Like Protein 5.
Immunogen Sequence	



Western blot analysis of lysates from THP-1 cells. primary antibody was diluted at 1:1000, 4°C over night

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081