

Anti-WFIKKN2 antibody (STJ194572)

STJ194572

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Wap-Kazal-Immunoglobulin-Kunitz And Ntr Domain-Containing Protein 2 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	124857
Gene Symbol	WFIKKN2
Uniprot ID	WFKN2_HUMAN
Immunogen	Synthesized peptide derived from human WFKN2
Region	
Specificity	WFIKKN2 polyclonal antibody (Wap-Kazal-Immunoglobulin-Kunitz And Ntr Domain-Containing Protein 2) binds to endogenous Wap-Kazal-Immunoglobulin-Kunitz And Ntr Domain-Containing Protein 2.
Immunogen Sequence	



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4°C over night

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081