

Anti-DEGS2 antibody (STJ194812)

STJ194812

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Sphingolipid Delta (4-Desaturase/C4-Monooxygenase Des2) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	123099
Gene Symbol	DEGS2
Uniprot ID	DEGS2_HUMAN
Immunogen	Synthesized peptide derived from human DEGS2
Immunogen Region	
Specificity	DEGS2 polyclonal antibody (Sphingolipid Delta (4-Desaturase/C4-Monooxygenase Des2) binds to endogenous Sphingolipid Delta (4-Desaturase/C4-Monooxygenase Des2).
Immunogen Sequence	



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4°C over night

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081