

Anti-NUDT6 antibody (STJ194964)

STJ194964

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Nucleoside Diphosphate-Linked Moiety X Motif 6 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11162
Gene Symbol	NUDT6
Uniprot ID	NUDT6_HUMAN
Immunogen	Synthesized peptide derived from human NUDT6
Immunogen Region	
Specificity	NUDT6 polyclonal antibody (Nucleoside Diphosphate-Linked Moiety X Motif 6) binds to endogenous Nucleoside Diphosphate-Linked Moiety X Motif 6.
Immunogen Sequence	



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°C over night