

Anti-FBXL19-AS1 antibody (STJ194981)

STJ194981

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Putative Uncharacterized Protein Fbxl19-As1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	FBXL19-AS1
Uniprot ID	FBAS1_HUMAN
Immunogen	Synthesized peptide derived from human FBAS1
Immunogen Region	
Specificity	FBXL19-AS1 polyclonal antibody (Putative Uncharacterized Protein Fbxl19-As1) binds to endogenous Putative Uncharacterized Protein Fbxl19-As1.
Immunogen Sequence	



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4°C over night