

Anti-RTN4IP1 antibody (STJ195183)

STJ195183

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Reticulon-4-Interacting Protein 1-Mitochondrial is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84816
Gene Symbol	RTN4IP1
Uniprot ID	RT4I1_HUMAN
Immunogen	Synthesized peptide derived from human RT4I1
Immunogen Region	
Specificity	RTN4IP1 polyclonal antibody (Reticulon-4-Interacting Protein 1-Mitochondrial) binds to endogenous Reticulon-4-Interacting Protein 1-Mitochondrial.
Immunogen Sequence	



Western blot analysis of lysates from HpeG2 cells. primary antibody was diluted at 1:1000, 4°C over night