

Anti-DENND2C antibody (STJ195645)

STJ195645

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Denn Domain-Containing Protein 2c is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	163259
Gene Symbol	DENND2C
Uniprot ID	DEN2C_HUMAN
Immunogen	Synthesized peptide derived from human DEN2C
Immunogen Region	
Specificity	DENND2C polyclonal antibody (Denn Domain-Containing Protein 2c) binds to endogenous Denn Domain-Containing Protein 2c.
Immunogen Sequence	



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4°C over night