

Anti-ANKRD24 antibody (STJ195770)

STJ195770

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ankyrin Repeat Domain-Containing Protein 24 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	170961
Gene Symbol	ANKRD24
Uniprot ID	ANR24_HUMAN
Immunogen	Synthesized peptide derived from human ANR24
Immunogen Region	
Specificity	ANKRD24 polyclonal antibody (Ankyrin Repeat Domain-Containing Protein 24) binds to endogenous Ankyrin Repeat Domain-Containing Protein 24.
Immunogen Sequence	



Western blot analysis of lysates from HpeG2 cells. primary antibody was diluted at 1:1000, 4°C over night