

Anti-ALOX15 antibody (STJ195947)

STJ195947

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Polyunsaturated Fatty Acid Lipoxygenase Alox15 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	246
Gene Symbol	ALOX15
Uniprot ID	LOX15_HUMAN
Immunogen	Synthesized peptide derived from human LOX15
Immunogen Region	
Specificity	ALOX15 polyclonal antibody (Polyunsaturated Fatty Acid Lipoxygenase Alox15) binds to endogenous Polyunsaturated Fatty Acid Lipoxygenase Alox15.
Immunogen Sequence	



Western blot analysis of lysates from THP-1 cells. primary antibody was diluted at 1:1000, 4°C over night

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081