

Anti-ANKRD27 antibody (STJ196146)

STJ196146

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ankyrin Repeat Domain-Containing Protein 27 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84079
Gene Symbol	ANKRD27
Uniprot ID	ANR27_HUMAN
Immunogen	Synthesized peptide derived from human ANR27
Immunogen Region	
Specificity	ANKRD27 polyclonal antibody (Ankyrin Repeat Domain-Containing Protein 27) binds to endogenous Ankyrin Repeat Domain-Containing Protein 27.
Immunogen Sequence	



Western blot analysis of lysates from HEK293 cells. primary antibody was diluted at 1:1000, 4°C over night