

Anti-BMI1 antibody (STJ196151)

STJ196151

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Polycomb Complex Protein Bmi-1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100532731/648
Gene Symbol	BMI1
Uniprot ID	BMI1_HUMAN
Immunogen	Recombinant Protein derived from human Bmi-1
Immunogen Region	
Specificity	BMI1 polyclonal antibody (Polycomb Complex Protein Bmi-1) binds to endogenous Polycomb Complex Protein Bmi-1.
Immunogen Sequence	