

Anti-PIK3R1 antibody (STJ196152)

STJ196152

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phosphatidylinositol 3-Kinase Regulatory Subunit Alpha is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5295
Gene Symbol	PIK3R1
Uniprot ID	P85A_HUMAN
Immunogen	Recombinant Protein derived from human PI3 Kinase P85
Immunogen Region	
Specificity	PIK3R1 polyclonal antibody (Phosphatidylinositol 3-Kinase Regulatory Subunit Alpha) binds to endogenous Phosphatidylinositol 3-Kinase Regulatory Subunit Alpha.
Immunogen Sequence	