

Anti-CALR antibody (STJ196154)

STJ196154

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Calreticulin is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	811
Gene Symbol	CALR
Uniprot ID	CALR_HUMAN
Immunogen	Recombinant Protein derived from human Calreticulin
Immunogen Region	
Specificity	CALR polyclonal antibody (Calreticulin) binds to endogenous Calreticulin.
Immunogen Sequence	