

Anti-KNG1 antibody (STJ196155)

STJ196155

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Kininogen-1 is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3827
Gene Symbol	KNG1
Uniprot ID	KNG1_HUMAN
Immunogen	Recombinant Protein derived from human Kininogen 1
Immunogen Region	
Specificity	KNG1 polyclonal antibody (Kininogen-1) binds to endogenous Kininogen-1.
Immunogen Sequence	