

Anti-TPO antibody (STJ196166)

STJ196166

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Thyroid Peroxidase is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7173
Gene Symbol	TPO
Uniprot ID	PERT_HUMAN
Immunogen	Synthesized peptide derived from human Thyroid Peroxidase (TPO)
Immunogen Region	
Specificity	TPO monoclonal antibody (Thyroid Peroxidase) binds to endogenous Thyroid Peroxidase.
Immunogen Sequence	