

Anti-Mast Cell Tryptase antibody [PT2025] (STJ196170)

STJ196170

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Mast cell tryptase is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PT2025
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthesized peptide derived from human Mast cell tryptase
Immunogen Region	
Specificity	Mast Cell Tryptase monoclonal antibody (Mast cell tryptase) binds to endogenous Mast cell tryptase.
Immunogen Sequence	