

Anti-ANKRD11 antibody (STJ196173)

STJ196173

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ankyrin Repeat Domain-Containing Protein 11 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29123
Gene Symbol	ANKRD11
Uniprot ID	ANR11_HUMAN
Immunogen	Synthesized peptide derived from human ANR11
Immunogen Region	
Specificity	ANKRD11 polyclonal antibody (Ankyrin Repeat Domain-Containing Protein 11) binds to endogenous Ankyrin Repeat Domain-Containing Protein 11.
Immunogen Sequence	