

Anti-LYRM4 antibody (STJ196186)

STJ196186

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Lyr Motif-Containing Protein 4 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57128
Gene Symbol	LYRM4
Uniprot ID	LYRM4_HUMAN
Immunogen	Synthesized peptide derived from human LYRM4
Immunogen Region	
Specificity	LYRM4 polyclonal antibody (Lyr Motif-Containing Protein 4) binds to endogenous Lyr Motif-Containing Protein 4.
Immunogen Sequence	