

Anti-SPINK4 antibody (STJ196189)

STJ196189

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine Protease Inhibitor Kazal-Type 4 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27290
Gene Symbol	SPINK4
Uniprot ID	ISK4_HUMAN
Immunogen	Synthesized peptide derived from human ISK4
Immunogen Region	
Specificity	SPINK4 polyclonal antibody (Serine Protease Inhibitor Kazal-Type 4) binds to endogenous Serine Protease Inhibitor Kazal-Type 4.
Immunogen Sequence	