

Anti-MORN2 antibody (STJ196191)

STJ196191

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Morn Repeat-Containing Protein 2 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	729967
Gene Symbol	MORN2
Uniprot ID	MORN2_HUMAN
Immunogen	Synthesized peptide derived from human MORN2
Immunogen Region	
Specificity	MORN2 polyclonal antibody (Morn Repeat-Containing Protein 2) binds to endogenous Morn Repeat-Containing Protein 2.
Immunogen Sequence	