

Anti-MCRIP1 antibody (STJ196193)

STJ196193

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Mapk-Regulated Corepressor-Interacting Protein 1 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	348262
Gene Symbol	MCRIP1
Uniprot ID	MCR11_HUMAN
Immunogen	Synthesized peptide derived from human F195B
Immunogen Region	
Specificity	MCRIP1 polyclonal antibody (Mapk-Regulated Corepressor-Interacting Protein 1) binds to endogenous Mapk-Regulated Corepressor-Interacting Protein 1.
Immunogen Sequence	