

Anti-FAM167A-AS1 antibody (STJ196207)

STJ196207

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Uncharacterized Protein Fam167a-As1 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	FAM167A-AS1
Uniprot ID	FAAS1_HUMAN
Immunogen	Synthesized peptide derived from human CH012
Immunogen Region	
Specificity	FAM167A-AS1 polyclonal antibody (Uncharacterized Protein Fam167a-As1) binds to endogenous Uncharacterized Protein Fam167a-As1.
Immunogen Sequence	