

## Anti-COA4 antibody (STJ196209)

STJ196209

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Cytochrome C Oxidase Assembly Factor 4 Homolog-Mitochondrial is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	IHC-P, IF-P
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	IHC-P 1:50-200
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">51287</a>
<b>Gene Symbol</b>	<a href="#">COA4</a>
<b>Uniprot ID</b>	<a href="#">COA4_HUMAN</a>
<b>Immunogen</b>	Synthesized peptide derived from human CHCH8
<b>Immunogen Region</b>	
<b>Specificity</b>	COA4 polyclonal antibody (Cytochrome C Oxidase Assembly Factor 4 Homolog-Mitochondrial) binds to endogenous Cytochrome C Oxidase Assembly Factor 4 Homolog-Mitochondrial.
<b>Immunogen Sequence</b>	