

Anti-EYS antibody (STJ196212)

STJ196212

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein Eyes Shut Homolog is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	346007
Gene Symbol	EYS
Uniprot ID	EYS_HUMAN
Immunogen	Synthesized peptide derived from human EYS
Immunogen Region	
Specificity	EYS polyclonal antibody (Protein Eyes Shut Homolog) binds to endogenous Protein Eyes Shut Homolog.
Immunogen Sequence	