

Anti-Phospho-TP53BP1-Thr543 antibody (STJ196235)

STJ196235

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Tp53-Binding Protein 1-Thr543 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7158
Gene Symbol	TP53BP1
Uniprot ID	TP53B_HUMAN
Immunogen	Synthesized phospho peptide around human 53BP1 (Thr543)
Immunogen Region	
Specificity	Phospho-TP53BP1-Thr543 polyclonal antibody (Tp53-Binding Protein 1) binds to endogenous Tp53-Binding Protein 1 only when phosphorylated at Thr543.
Immunogen Sequence	