

Anti-Phospho-ALK-Tyr1278 antibody (STJ196237)

STJ196237

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Alk Tyrosine Kinase Receptor-Tyr1278 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	238
Gene Symbol	ALK
Uniprot ID	ALK_HUMAN
Immunogen	Synthesized phospho peptide around human ALK (Tyr1278)
Immunogen Region	
Specificity	Phospho-ALK-Tyr1278 polyclonal antibody (Alk Tyrosine Kinase Receptor) binds to endogenous Alk Tyrosine Kinase Receptor only when phosphorylated at Tyr1278.
Immunogen Sequence	