

Anti-Phospho-ARHGAP42-Tyr376 antibody (STJ196241)

STJ196241

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Rho Gtpase-Activating Protein 42-Tyr376 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	143872
Gene Symbol	ARHGAP42
Uniprot ID	RHG42_HUMAN
Immunogen	Synthesized phosho peptide around human ARHGAP42 (Tyr376)
Immunogen Region	
Specificity	Phospho-ARHGAP42-Tyr376 polyclonal antibody (Rho Gtpase-Activating Protein 42) binds to endogenous Rho Gtpase-Activating Protein 42 only when phosphorylated at Tyr376.
Immunogen Sequence	