

Anti-Phospho-ATG13-Ser355 antibody (STJ196246)

STJ196246

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Autophagy-Related Protein 13-Ser355 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9776
Gene Symbol	ATG13
Uniprot ID	ATG13_HUMAN
Immunogen	Synthesized phosho peptide around human Atg13 (Ser355)
Immunogen Region	
Specificity	Phospho-ATG13-Ser355 polyclonal antibody (Autophagy-Related Protein 13) binds to endogenous Autophagy-Related Protein 13 only when phosphorylated at Ser355.
Immunogen Sequence	