

Anti-Phospho-BCL9L-Ser915 antibody (STJ196251)

STJ196251

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-B-Cell Cll/Lymphoma 9-Like Protein-Ser915 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	283149
Gene Symbol	BCL9L
Uniprot ID	BCL9L_HUMAN
Immunogen	Synthesized phospho peptide around human BCL9L (Ser915)
Immunogen Region	
Specificity	Phospho-BCL9L-Ser915 polyclonal antibody (B-Cell Cll/Lymphoma 9-Like Protein) binds to endogenous B-Cell Cll/Lymphoma 9-Like Protein only when phosphorylated at Ser915.
Immunogen Sequence	