

Anti-Phospho-CBL-Tyr731 antibody (STJ196267)

STJ196267

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-E3 Ubiquitin-Protein Ligase Cbl-Tyr731 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	867
Gene Symbol	CBL
Uniprot ID	CBL_HUMAN
Immunogen	Synthesized phospho peptide around human c-Cbl (Tyr731)
Immunogen Region	
Specificity	Phospho-CBL-Tyr731 polyclonal antibody (E3 Ubiquitin-Protein Ligase Cbl) binds to endogenous E3 Ubiquitin-Protein Ligase Cbl only when phosphorylated at Tyr731.
Immunogen Sequence	